


Cadenza CW-C18

Technical

Basic Interaction of Cadenza CW-C18



Ref) N.Tanaka et.al., J. Chromatg. Sci., 27, 721 (1989)

Cadenza CW-C18 is a high-resolution and high-throughput ODS column which consists of 30nm pore and 3um particle size. Due to the larger pore size, there are several key advantages:

- (1) Better end-capping efficiency provides improved peak shape
- (2) High-speed analysis under the same conditions as conventional ODS column
- (3) Improvement of peak response for relatively large molecules

It may provide different selectivity from other ODS phases due to the difference of molecular interaction as follows:

**[Hydrophobicity]**

Hydrophobicity is a fundamental property of ODS phase. CW-C18 results in almost the same value which means the same ODS ligand density as CD-C18 (same phase structure).

**[Electrostatic Interaction]**

This interaction is related to the existence of siloxane and silanol on the surface of silica material. CW-C18 showed a minimum value and was similar to CD-C18. These low values come from a unique and effective end-capping technology, "Polymeric-end-capping."

**[Steric Selectivity]**

High ligand density ODS phase has the ability to recognize steric structure of similar compounds. CW-C18 showed a large value similar to CD-C18. CW-C18 is useful for high-throughput analysis with excellent molecular recognition.