


Intrada Organic Acid

150 x 2 mm

Application

LC-MS: Oxalic acid (sample solvent effect)

シュウ酸のLC-MS分析 (試料溶媒の影響)



Intrada Organic Acid, 150 x 2 mm

A: acetonitrile / water / HCOOH = 10 / 90 / 0.1

B: acetonitrile / 100mM HCOONH₄ = 10 / 90

0%B (0-1min), 0-100%B (1-7min), 100%B (7-10min)

0.2mL/min (12MPa), 37°C, 5µL

MS (ESI, SIM), [M-H]⁻ m/z = 89.0