


Nardis ND-NX

250 x 4.6 mm

HPLC COLUMN APPLICATION

Group separation of monosaccharides

単糖類のグループ分離



Nardis ND-NX, 250 x 4.6 mm

A: acetonitrile / 100mM ammonium acetate = 95 / 5

B: acetonitrile / water / acetic acid = 10 / 90 / 0.8

0 -100%B (0 -20 min), 100%B (20 -30 min),

0.8 mL/min (5MPa), 60°C

ELSD (spray chamber 70°C, drift chamber 100°C)

10 µL (0.5-1.5mg/mL, water)