
Mintakt Rechnical Information No.TI371E



Unison UK-C18, 150 x 4.6 mm

A: 50mM sodium phosphate buffer (pH 3.2) / acetonitrile = 91 / 9 (v/v) B: acetonitrile 0%B (0-18 min), 0-30%B (18-28 min), 30%B (28-32 min) 0.75 mL/min, 40 deg.C, Ex.470 nm / Em. 530 nm

Pre-column derivatization: 105 ul of 100mM borate buffer (pH 9.0) and 30 uL of 40mM NBD-F in acetonitrile were added to the 30 uL of sample. Then, the reaction was incubated at 40 deg.C for 3 min. To stop the reaction, 435 uL of 0.5% acetic acid (v/v) was added. A10 uL aliquot of sample was injected onto the HPLC.

Journal of Chromatography A, 1066 (2005) 41-45

Courtesy of Dr. TSUNODA, University of Tokyo, JAPAN