


Intrada Organic Acid

150 x 2 mm

HPLC COLUMN APPLICATION

LC-MS: Adenosine phosphates

アデノシンリン酸の LCMS 分析



Intrada Organic Acid, 150 x 2 mm

A: acetonitrile / 100mM ammonium acetate = 95 / 5

B: acetonitrile / 100mM ammonium acetate = 10 / 90

0-60%B (0-11min), 100%B (11-16min)

0.2mL/min (5MPa), 60°C, MS (ESI, SIM, Negative)

5µL (50µg/mL, ACN / water / HCOOH = 50 / 50 / 1)