


Unison UK-C1  
 Unison UK-Phenyl Cadenza CX-C18  
 Unison UK-C8 Cadenza CW-C18  
 Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm


HPLC COLUMN APPLICATION

## (2-Hydroxy-5-methylphenyl)triphenylphosphonium bromide (Condensation organocatalysts)

(2-ヒドロキシ-5-メチルフェニル)トリフェニルホスホニウムブロミド  
 (縮合有機触媒)



(2-hydroxy-5-methylphenyl)  
 triphenylphosphonium bromide



100 x 4.6 mm

A: water / formic acid = 100 / 0.5

B: water / acetonitrile / formic acid = 10 / 90 / 0.5

0%B (0-1min), 0-100%B (1-15min)

1.0mL/min (10MPa), 37°C, 260nm

10µL (0.5mg/mL, water / acetonitrile / formic acid = 90 / 10 / 1)