
Mintakt Rechnical Information No.TI166E

Unison UK-Silica150 x 4.6 mmTechnicalHigh-efficiency separation by a 3um silica particle



The Unison series of columns provides authoritative results in normal phase separation using its silica stationary phase. This high level of performance is in addition to its extraordinary resolution, higher than conventional 5um columns.

The example above shows a normal phase separation of the tocopherol isomers. Beta- and gammatocopherol are position isomers, compounds characterized as difficult to separate. Unison UK-Silica 150 x 4.6mm, however, provides excellent separation in under 10 minutes. By using the 3um silica Unison column, you can improve your efficiency and analysis speed compared with conventional 5um, 25cm columns.

The Unison UK-Silica offers optimal solutions for specialty research areas including natural chemistry and organic combinatorial chemistry.

USA, Canada, and Mexico : Silvertone Sciences (USA) - info@silvertonesciences.com Other Countries : Imtakt Corp.(JAPAN) - info@imtakt.com