

Presto FF-C18

250 x 4.6 mm

Application

**Recombinant and biosynthetic human serum albumin (HSA)**

遺伝子組換えおよび生合成によるヒト血清アルブミン(HSA)

Presto FF-C18, 250 x 4.6 mm

A: water / trifluoroacetic acid = 100 / 0.1

B: acetonitrile / trifluoroacetic acid = 100 / 0.1

35-41 %B (0-120 min)

0.4 mL/min (17 MPa), 80°C, 220 nm

20 µL (10mg/mL, 0.1% TFA)

*recombinant HSA**biosynthetic HSA*