


ODS Surface Structures


Imtakt's 3µm particle technology

Unison UK-C18




Optimized ODS coverage
Fully polymeric end-capping
100% aqueous elution
100% organic solvent elution

Cadenza CD-C18
Cadenza CW-C18




Fully polymeric end-capping
Superior molecular recognition

Cadenza CL-C18



controlled residual silanols
100% aqueous elution
Different selectivity from CD-C18

Cadenza HS-C18



Hybrid phase with
ODS and polar ligands
End-capping
Ultra low-hydrophobicity
Different selectivity