

Intrada SEC

250 x 4.6 mm

Application

SEC: Ribonuclease A (RNase A), from bovine pancreas

リボヌクレアーゼAのSEC分離

Ribonuclease A (RNase A)
from bovine pancreas
MW 13kDa

ELSD

0 5 10 15 min

UV 280nm

0 5 10 15 min

Intrada SEC, 250 x 4.6 mm
100mM ammonium formate
0.3 mL/min (5MPa), 37°C
10 µL (1mg/mL, 10mM ammonium acetate)