




Unison Surface Structures


Imtakt's 3um particle technology

Unison UK-C18




*Optimized ODS coverage
Fully polymeric end-capping
100% aqueous elution
100% organic solvent elution*

Unison UK-C8




*Optimized C8 coverage
Fully polymeric end-capping
100% aqueous elution
100% organic solvent elution*

Unison UK-Phenyl



*Optimized Phenyl coverage
Fully polymeric end-capping
100% aqueous elution
100% organic solvent elution*

Unison UK-Silica



*Siloxane / Silanol ligand
100% aqueous acidic elution
100% organic solvent elution*