
Unison UK-Phenyl Unison UK-C1
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm HPLC COLUMN APPLICATION

Pongamol (Natural ultraviolet absorber)

ポンガモール（天然紫外線吸収剤）



pongamol

100 x 4.6 mm
A: water / acetonitrile / acetic acid = 90 / 10 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
45-100%B (0-15min)
1.0mL/min (8-12MPa), 37°C, 350nm
10µL (0.025mg/mL, water / acetonitrile / acetic acid = 20 / 80 / 1)