
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm **HPLC COLUMN APPLICATION**

Glabridin (Isoflavonoids)

グラブリジン（イソフラボノイド類）



glabridin

100 x 4.6 mm
A: water / acetonitrile / acetic acid = 90 / 10 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
35-100%B (0-15min)
1.0mL/min (8-12MPa), 37°C, 260nm
10µL (0.2mg/mL, water / acetonitrile / acetic acid = 70 / 30 / 1)