
Cadenza *CX-C18* 150 x 3 mm **HPLC COLUMN APPLICATION**

(2S)-N-[(1S)-1-(Hydroxydiphenylmethyl)-3-methylbutyl] -2-pyrrolidinecarboxamide (Proline-organocatalyst)

(2S)-N-[(1S)-1-(ヒドロキシジフェニルメチル)-3-メチルブチル] -2-ピロリジンカルボキサミド (プロリン型有機触媒)



(2S)-N-[(1S)-1-(hydroxydiphenylmethyl) -3-methylbutyl]-2-pyrrolidinecarboxamide

150 x 3 mm
A: water / acetonitrile / formic acid = 95 / 5 / 0.1
B: water / acetonitrile / formic acid = 50 / 50 / 0.1
0-75%B (0-20min)
0.4mL/min (11MPa), 37°C, 260nm
2µL (2.0mg/mL, water / formic acid = 100 / 1)