
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm HPLC COLUMN APPLICATION

Galangin (Flavonols)

ガランギン（フラボノール類）



galangin

100 x 4.6 mm
A: water / acetonitrile / acetic acid = 90 / 10 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
25-100%B (0-15min)
1.0mL/min (9-16MPa), 37°C, 360nm
10µL (0.05mg/mL, water / acetonitrile / acetic acid = 30 / 70 / 1)