
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm **HPLC COLUMN APPLICATION**

Morin (Flavonols)

モリン（フラボノール類）



100 x 4.6 mm
A: water / acetonitrile / formic acid = 90 / 10 / 0.5
B: water / acetonitrile / formic acid = 10 / 90 / 0.5
0-65%B (0-15min)
1.0mL/min (10-13MPa), 37°C, 260nm
10µL (0.05mg/mL, water / acetonitrile / formic acid = 80 / 20 / 1)