
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18


100 x 4.6 mm **HPLC COLUMN APPLICATION**

Morin (Flavonols)

モリン（フラボノール類）



morin



100 x 4.6 mm
A: 100mM ammonium formate / acetonitrile = 95 / 5
B: 100mM ammonium formate / acetonitrile = 10 / 90
0%B (0-1min), 0-50%B (1-15min)
1.0mL/min (10-14MPa), 37°C, 260nm
10µL (0.1mg/mL, 100mM ammonium formate / acetonitrile = 70 / 30)