
Cadenza CL-C18
Cadenza CX-C18
Unison UK-C18
Cadenza CD-C18


100 x 4.6 mm HPLC COLUMN APPLICATION

Phenothrin (anthelminthic agent)

フェノトリン（駆虫薬）



phenothrin



100 x 4.6 mm
A: water / acetonitrile / acetic acid = 50 / 50 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
25-100%B (0-15min), 100%B (15-20min)
1.0mL/min (8-12MPa), 37°C, 260nm
10µL (1.5mg/mL, water / acetonitrile / acetic acid = 40 / 60 / 1)