
Unison UK-Amino
Nardis ND-NX

100 x 4.6 mm HPLC COLUMN APPLICATION

(S)-5-(Pyrrolidin-2-yl)-1H-tetrazole (Proline-type organocatalyst)

(S)-5-(ピロリジン-2-イル)-1H-テトラゾール (プロリン型有機触媒)



(S)-5-(pyrrolidin-2-yl)-1H-tetrazole

0-100%B (0-15min)
1.0mL/min (5-6MPa), 37°C, ELSD (spray chamber 30°C, drift tube 60°C)
10µL (2.0mg/mL, 100mM ammonium formate / acetonitrile = 10 / 90)