
Cadenza CD-C18
Cadenza CW-C18
Unison UK-C18

100 x 3 mm

HPLC COLUMN APPLICATION

# LC-MS: Chloroquine analogues (Antimalarial)

クロロキン類縁体（マラリア薬）のLC-MS分析




*100 x 3 mm*
*A: 100mM ammonium acetate / acetonitrile = 95 / 5*
*B: 100mM ammonium acetate / acetonitrile = 5 / 95*
*0-60%B (0-15min)*
*0.4mL/min (8-10MPa), 37°C, Single Quad. MS (ESI)*
*5µL (2µg/mL, 100mM ammonium acetate)*