
Scherzo SM-C18
Nardis ND-RX 100 x 3 mm HPLC COLUMN APPLICATION

LC-MS: Chloroquine analogues (Antimalarial)

クロロキン類縁体（マラリア薬）の LC-MS分析



100 x 3 mm
A: 10mM ammonium formate / acetonitrile = 95 / 5
B: water / acetonitrile / formic acid = 50 / 50 / 0.5
0-80%B (0-15min)
0.4mL/min (7-8MPa), 37°C, Single Quad. MS (ESI)
5µL (1µg/mL, 100mM ammonium formate = 100 / 1)