
Cadenza CD-C18
Cadenza CW-C18
Unison UK-C18

100 x 3 mm


HPLC COLUMN APPLICATION

LC-MS: Oseltamivir phosphate (Anti-influenza virus agent)

オセルタミビル リン酸塩（抗インフルエンザウイルス薬）の LC-MS分析



oseltamivir phosphate



100 x 3 mm
A: 100mM ammonium formate / acetonitrile = 90 / 10
B: 100mM ammonium formate / acetonitrile = 5 / 95
0-100%B (0-10min)
0.4mL/min (8-10MPa), 37°C, Single Quad. MS (ESI)
5µL (1µg/mL, 100mM ammonium formate)