
Unison UK-Silica 150 x 3 mm
Unison UK-C18 100 x 3 mm

HPLC COLUMN APPLICATION


LC-MS: Oseltamivir phosphate (Anti-influenza virus agent)

オセルタミビル リン酸塩（抗インフルエンザウイルス薬）の LC-MS分析



oseltamivir phosphate

Unison UK-Silica
150 x 3 mm
A: ACN / HCOOH = 100 / 0.1
B: ACN / water / HCOOH = 20 / 80 / 0.5
0-100%B (0-15min)
0.4mL/min (3MPa), 37°C,
Single Quad. MS (ESI)
5µL (1µg/mL,
ACN /water /HCOOH = 50 /50 / 1)



Unison UK-C18
100 x 3 mm
A: 100mM HCOONH4 / ACN = 90 / 10
B: 100mM HCOONH4 / ACN = 5 / 95
0-100%B (0-10min)
0.4mL/min (8MPa), 37°C,
Single Quad. MS (ESI)
5µL (1µg/mL, HCOONH4)