
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm **HPLC COLUMN APPLICATION**

Phthalimide (Non-heterocycle building block)

フタルイミド　(非ヘテロ環ビルディングブロック)



phthalimide

100 x 4.6 mm
A: water / acetonitrile / acetic acid = 95 / 5 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
0%B (0-1min), 0-55%B (1-15min)
1.0mL/min (11-13MPa), 37°C, 260nm
10µL (0.5mg/mL, water / DMSO / acetic acid = 40 / 60 / 1)