
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm **HPLC COLUMN APPLICATION**

S-(+)-Methoprene (Pesticide)

S-(+)-メトプレン (農薬)



S-(+)-methoprene

100 x 4.6 mm
A: water / acetonitrile / acetic acid = 50 / 50 / 0.5
B: acetonitrile / acetic acid = 100 / 0.5
10-100%B (0-15min)
1.0mL/min (12-16MPa), 37°C, 260nm
10µL (3.0mg/mL, water / acetonitrile / acetic acid = 10 / 90 / 1)