
Unison UK-C1
Unison UK-Phenyl Cadenza *CX-C18*
Unison UK-C8 Cadenza CW-C18
Unison UK-C18 Cadenza CD-C18


100 x 4.6 mm **HPLC COLUMN APPLICATION**

Icariin (Flavonoid glycoside)

イカリイン (フラボノイド配糖体)



icariin



100 x 4.6 mm
A: water / acetonitrile / acetic acid = 90 / 10 / 0.5
B: water / acetonitrile / acetic acid = 10 / 90 / 0.5
0-75%B (0-15min)
1.0mL/min (12-14MPa), 37°C, 260nm
10µL (0.075mg/mL, water / acetonitrile / acetic acid = 40 / 60 / 1)